

# CHEMISTRY 235 EXPERIMENT 3

## QUALITATIVE ANALYSIS

### METHOD AND RESULTS

You are provided with two samples labelled A and B. Each one contains one of the ions  $\text{Al}^{3+}$ , and  $\text{Fe}^{2+}$ . Perform the following tests so as to identify the ion present in each substance. Write net ionic equations wherever appropriate.

#### COMPOUND A

<i>TESTS</i>	<i>OBSERVATIONS</i>	<i>DEDUCTIONS</i>
Use a freshly prepared soln. of A in water for the tests below. Use a fresh portion for each test unless otherwise instructed.		
a) Add $\text{NH}_3(\text{aq})$ until in excess.		
b) Add $\text{NaOH}(\text{aq})$ until in excess. Boil the result.		
c) Add dil. $\text{HCl}$ then $\text{H}_2\text{S}$		
d) Add $\text{CH}_3\text{COONa}(\text{aq})$ then $\text{H}_2\text{S}(\text{aq})$		
e) Add $\text{K}_4\text{Fe}(\text{CN})_6(\text{aq})$ . Divide the result into 2 parts i) Add $\text{H}_2\text{O}_2(\text{aq})$ to 1 part. ii) Leave 2nd part to stand for 5 minutes.		
f) Add $\text{K}_3\text{Fe}(\text{CN})_6(\text{aq})$		

g) Add $\text{NH}_4\text{SCN}(\text{aq})$		
h) Add butanedionedioxime $(\text{aq})^1$ + ammonia soln.		
Add $1 \text{ cm}^3$ of conc. $\text{HNO}_3$ to $5 \text{ cm}^3$ of $\text{A}(\text{aq})$ . Bring to boil. Use this soln. for the following tests.		
a) Add $\text{NaOH}(\text{aq})$ until in excess.		
b) Add $\text{NH}_3$ soln. until in excess.		
c) Add $\text{H}_2\text{S}$ soln.		
d) Add $\text{K}_4\text{Fe}(\text{CN})_6(\text{aq})$ dropwise.		
e) Add $\text{K}_3\text{Fe}(\text{CN})_6(\text{aq})$ dropwise		
f) Add $\text{NH}_4\text{SCN}(\text{aq})$		
g) Add butanedionedioxime $(\text{aq})$		

---

<sup>1</sup> Also known as dimethylglyoxime, DMG.

## COMPOUND B

<i>TESTS</i>	<i>OBSERVATION</i>	<i>DEDUCTION</i>
Use an aqueous solution of B for the following tests. Use a fresh portion for each test unless otherwise instructed.		
a) Add NaOH(aq) until in excess. Divide the result into 2 parts. i) Add NH <sub>4</sub> Cl(aq) to 1st part. ii) Add dil. HCl to 2nd part.		
b) Add NH <sub>4</sub> Cl(s) followed by NH <sub>3</sub> (aq) until in excess.		
c) Add dil. HCl then H <sub>2</sub> S(aq)		
d) Add H <sub>2</sub> S(aq)		
e) Add Na <sub>2</sub> CO <sub>3</sub> (aq)		
f) Add 1 or 2 drops of litmus soln., then dil. HCl until acid (if the soln. is not already acid), then add NH <sub>3</sub> (aq) until alkaline. Allow ppt. to settle and observe it.		