

CHEMISTRY 235 EXPERIMENT 6

QUALITATIVE ANALYSIS

METHOD AND RESULTS

You are provided with three salts containing the cations Zn^{2+} , Mg^{2+} , NH_4^+ . Identify the cation present in each salt by performing the following tests. *Make sure you carry out confirmatory tests on any gases evolved.* Write your deductions and equations (where appropriate) in the deduction columns.

COMPOUND A

TEST	OBSERVATIONS	DEDUCTIONS
a) Perform a flame test on A. Use a fresh portion for each test unless otherwise instructed.		
Prepare an aqueous solution of A and use it for the following tests.		
a) Add $\text{NaOH}(\text{aq})$ until in excess.		
b) Add $\text{NH}_3(\text{aq})$ until in excess.		
c) Add $\text{Na}_2\text{CO}_3(\text{aq})$ then dil. HCl or dil. HNO_3 .		
d) i) Add $\text{NH}_4\text{Cl}(\text{aq})$ then $\text{NH}_3(\text{aq})$ then $\text{Na}_3\text{PO}_4(\text{aq})$. ii) Add dil. HCl or HNO_3 to the ppt. from (i).		
e) Add $\text{Na}_2(\text{COO})_2(\text{aq})$.		

TEST	OBSERVATIONS	DEDUCTIONS
f) Add $K_2CrO_4(aq)$.		
g) Add H_2S soln.		

COMPOUND B

TEST	OBSERVATIONS	DEDUCTIONS
Make a solution of B and use it for the following tests. Use a fresh portion for each test unless otherwise instructed.		
a) Perform a flame test on B.		
b) Add $NaOH(aq)$ until in excess.		
c) Add $NH_3(aq)$ until in excess.		
d) Add H_2S soln. then dil. HCl .		
e) Add $K_4Fe(CN)_6(aq)$ then $NH_3(aq)$.		
f) Add $K_3Fe(CN)_6(aq)$ then $NH_3(aq)$.		

COMPOUND C

TEST	OBSERVATIONS	DEDUCTIONS
a) Heat a small portion of C in a test tube.		
b) Add NaOH(aq) to a soln. of C and warm.		